

Product Datasheet

Murine IL-12 ELISA Kit 660.040

Product Name	Murine IL-12 ELISA Kit
Catalog Number	660.040
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-2000 pg/ml
Sensitivity	5.9 pg/ml
Tested Applications	ELISA
Alternative Names	IL-12
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant murine IL-12
Sample Types	Serum, Plasma, Cell culture supernatant