

## Product Datasheet

### Goat Immunoglobulin G (IgG) H + L Antibody A151

**Product Name** Goat Immunoglobulin G (IgG) H + L Antibody

**Catalog Number** A151

**Species/Host** Rabbit

**Tested Applications** ELISA, WB, IEP

**Immunogen** Purified Goat IgG

**Brand** EastCoast Bio

**Form/Appearance** Affinity Purified Polyclonal Antibody

**Storage** 2-8°C. Do not freeze.

**Clonality** Polyclonal

**Purity** >95%

**Buffer** 10 mM Phosphate Buffered Saline, pH 7.2

**Preservatives** 0.1% Sodium Azide