

## Product Datasheet

### Parvovirus, Canine Antibody A312

<b>Product Name</b>	Parvovirus, Canine Antibody
<b>Catalog Number</b>	A312
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, WB, IHC, IFA
<b>Immunogen</b>	Canine Parvovirus
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Reactivity with canine parvovirus, feline panleukopenia virus and mink parvovirus. No reactivity with canine adenovirus type 2, canine corona virus, canine distemper virus, canine parainfluenza virus, feline leukemia virus and feline immunodeficiency virus.
<b>Isotype</b>	IgG2a, kappa
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Cross Reactivity</b>	Reactivity with Pig and Mink Parvoviruses. No reactivity with Canine Adenovirus type 2, Canine Coronavirus, Canine Distemper Virus, Canine Parainfluenza Virus, Feline Leukemia Virus and Feline Immunodeficiency Virus
<b>Specificity</b>	Canine Parvovirus and Feline Panleucopenia Virus
<b>Matched Pairs</b>	A312 (self-pairs)
<b>Preservatives</b>	0.05% Sodium Azide