

Product Datasheet

Methamphetamine Antibody C03-99-31M-P

Product Name Methamphetamine Antibody

Catalog Number C03-99-31M-P

Species/Host Mouse

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Protein A Purified Monoclonal Antibody

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG2a, kappa

Clonality Monoclonal

Buffer 10 mM Phosphate Buffered Saline, pH 7.2

Specificity Methamphetamine

Matched Pairs P81-99-31A (capture) or P91-99-31H (detection)

Preservatives 0.05% Sodium Azide