

## Product Datasheet

### Amphetamine Antibody C03-99-72M-P

<b>Product Name</b>	Amphetamine Antibody
<b>Catalog Number</b>	C03-99-72M-P
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2b, kappa
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Specificity</b>	Amphetamine
<b>Matched Pairs</b>	P81-99-71A (capture) or P91-99-71H (detection)
<b>Preservatives</b>	0.05% Sodium Azide