

## Product Datasheet

### Amphetamine Antibody C03-99-72M-P

**Product Name** Amphetamine Antibody

**Catalog Number** C03-99-72M-P

**Species/Host** Mouse

**Tested Applications** LF, ELISA

**Brand** EastCoast Bio

**Form/Appearance** Protein A Purified Monoclonal Antibody

**Storage** Short Term ( $\leq 2$  weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Isotype** IgG2b, kappa

**Clonality** Monoclonal

**Buffer** 10 mM Phosphate Buffered Saline, pH 7.2

**Specificity** Amphetamine

**Matched Pairs** P81-99-71A (capture) or P91-99-71H (detection)

**Preservatives** 0.05% Sodium Azide