

Product Datasheet

Murine IL-22 ELISA Kit 660.110

| | |
|----------------------------|--|
| Product Name | Murine IL-22 ELISA Kit |
| Catalog Number | 660.110 |
| Description | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure. |
| Range | 7.8-500pg/ml |
| Sensitivity | 1.6pg/ml |
| Tested Applications | ELISA |
| Alternative Names | NULL |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Specificity | Recognizes both natural and recombinant murine IL-22 |
| Sample Types | Cell culture supernatant and serum tested, other body fluids suitable |