

Product Datasheet

Cytomegalovirus (CMV) Antigen, Native CV216

Product Name	Cytomegalovirus (CMV) Antigen, Native
Catalog Number	CV216
Tested Applications	ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Liquid
Storage	≤-70°C
Note	Confirmed negative for HIV and Hepatitis B Virus.
Purity	Lectin affinity purified glycoproteins
Source	NHDF Cells/Strain AD169
Buffer	50 mM Tris, 1 M NaCl, 1 mM MgCl2, 0.1% Triton X-100, 0.2 M Methylmannopyrannoside, pH 7.5