

Product Datasheet

Epstein-Barr Virus (EBV-gp110) Antibody EA068

Product Name	Epstein-Barr Virus (EBV-gp110) Antibody
Catalog Number	EA068
Species/Host	Mouse
Tested Applications	ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	gp110