

---

## Product Datasheet

### Goat anti-Mouse Immunoglobulin G (IgG) H + L Antibody F101

**Product Name** Goat anti-Mouse Immunoglobulin G (IgG) H + L Antibody

**Catalog Number** F101

**Species/Host** Goat

**Tested Applications** ELISA, WB

**Brand** EastCoast Bio

**Form/Appearance** Affinity Purified Polyclonal Antibody

**Storage** 2-8°C

**Clonality** Polyclonal

**Purity** ≥90%

**Buffer** 10 mM Phosphate Buffered Saline, pH 7.2

**Specificity** All IgG subclasses and all mouse Ig

**Preservatives** 0.09% Sodium Azide