
Product Datasheet

Goat anti-Mouse Immunoglobulin G (IgG) H + L Antibody F101

Product Name Goat anti-Mouse Immunoglobulin G (IgG) H + L Antibody

Catalog Number F101

Species/Host Goat

Tested Applications ELISA, WB

Brand EastCoast Bio

Form/Appearance Affinity Purified Polyclonal Antibody

Storage 2-8°C

Clonality Polyclonal

Purity ≥90%

Buffer 10 mM Phosphate Buffered Saline, pH 7.2

Specificity All IgG subclasses and all mouse Ig

Preservatives 0.09% Sodium Azide