

Product Datasheet

Thyroxine (T4) Antibody FA303

Product Name	Thyroxine (T4) Antibody
Catalog Number	FA303
Species/Host	in vitro cell culture
Tested Applications	ELISA
Immunogen	T4-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Mouse Monoclonal Antibody
Storage	2-8°C
Isotype	IgG 2b, kappa
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	T4
Preservatives	0.1% Sodium Azide