

## Product Datasheet

### Murine MCP-1 ELISA Kit 660.130

<b>Product Name</b>	Murine MCP-1 ELISA Kit
<b>Catalog Number</b>	660.130
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	15.6-1000pg/ml
<b>Sensitivity</b>	2.2pg/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	CCL2, SCYA-2, JE
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Specificity</b>	Recognizes both natural and recombinant murine MCP-1
<b>Sample Types</b>	Cell culture supernatant and serum tested, other body fluids suitable