

## Product Datasheet

### Heartworm, Canine Antibody FA306

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Heartworm, Canine Antibody               |
| <b>Catalog Number</b>      | FA306                                    |
| <b>Species/Host</b>        | Mouse                                    |
| <b>Tested Applications</b> | ELISA                                    |
| <b>Immunogen</b>           | R. G. antigen                            |
| <b>Brand</b>               | EastCoast Bio                            |
| <b>Form/Appearance</b>     | Purified Monoclonal Antibody             |
| <b>Storage</b>             | 2-8°C. Do not freeze.                    |
| <b>Isotype</b>             | IgG1                                     |
| <b>Clonality</b>           | Monoclonal                               |
| <b>Purity</b>              | >95%                                     |
| <b>Buffer</b>              | 0.05 M Phosphate Buffered Saline, pH 7.4 |
| <b>Specificity</b>         | Dirofilaria immitis                      |
| <b>Matched Pairs</b>       | JP051                                    |
| <b>Preservatives</b>       | 0.1% Sodium Azide                        |