

## Product Datasheet

### Human Immunoglobulin E (IgE) Antibody FA317

**Product Name** Human Immunoglobulin E (IgE) Antibody

**Catalog Number** FA317

**Species/Host** Mouse

**Tested Applications** ELISA

**Immunogen** Human IgE

**Brand** EastCoast Bio

**Form/Appearance** Purified Monoclonal Antibody

**Storage** 2-8°C. Do not freeze.

**Isotype** IgG2b

**Clonality** Monoclonal

**Purity** >95%

**Buffer** 0.05 M Phosphate Buffered Saline, pH 7.4

**Specificity** Human IgE

**Matched Pairs** FA318 (capture)

**Preservatives** 0.1% Sodium Azide