

Product Datasheet

Adenovirus Antibody FG300

Product Name	Adenovirus Antibody
Catalog Number	FG300
Species/Host	Mouse
Tested Applications	LF, ELISA, WB, IFA
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.0-7.5
Specificity	Hexon Protein
Matched Pairs	HM1178 (detection)