

Product Datasheet

Human Chorionic Gonadotropin α (HCG α) Antibody FK306

Product Name	Human Chorionic Gonadotropin α (HCG α) Antibody
Catalog Number	FK306
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Purified alpha HCG
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (\leq 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
Specificity	HCG α
Preservatives	0.05% Sodium Azide