

## Product Datasheet

### Human Chorionic Gonadotropin $\alpha$ (HCG $\alpha$ ) Antibody FK306

<b>Product Name</b>	Human Chorionic Gonadotropin $\alpha$ (HCG $\alpha$ ) Antibody
<b>Catalog Number</b>	FK306
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Purified alpha HCG
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq$ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
<b>Specificity</b>	HCG $\alpha$
<b>Preservatives</b>	0.05% Sodium Azide