

## Product Datasheet

### Influenza A Virus Antibody FT351

<b>Product Name</b>	Influenza A Virus Antibody
<b>Catalog Number</b>	FT351
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, WB, IHC
<b>Immunogen</b>	Influenza A Virus
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq$ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.0
<b>Cross Reactivity</b>	No reactivity with Influenza B, Adenovirus, and RSV
<b>Specificity</b>	Nucleoprotein from all strains
<b>Preservatives</b>	0.02% Sodium Azide