

Product Datasheet

Cancer Antigen 15-3 (CA 15-3) MUC 1 Antibody FT372

Product Name	Cancer Antigen 15-3 (CA 15-3) MUC 1 Antibody
Catalog Number	FT372
Species/Host	Mouse
Tested Applications	ELISA, IHC
Immunogen	Human milk fat globules
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (\leq 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffer, pH 7.0
Specificity	MUC1
Preservatives	0.02% Sodium Azide