

Product Datasheet

Cortisol-3-BSA Antigen FV501

Product Name Cortisol-3-BSA Antigen

Catalog Number FV501

Tested Applications LF

Brand EastCoast Bio

Form/Appearance Liquid; Antigen-BSA conjugate

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Purity $\geq 90\%$ by SDS-PAGE

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Preservatives 0.02% Sodium Azide