

## Product Datasheet

### Cortisol-3-BSA Antigen FV501

<b>Product Name</b>	Cortisol-3-BSA Antigen
<b>Catalog Number</b>	FV501
<b>Tested Applications</b>	LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Antigen-BSA conjugate
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	$\geq 90\%$ by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Preservatives</b>	0.02% Sodium Azide