

Product Datasheet

Cortisol-3-BSA Antigen FV501

Product Name	Cortisol-3-BSA Antigen
Catalog Number	FV501
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	$\geq 90\%$ by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservatives	0.02% Sodium Azide