

---

## Product Datasheet

### Testosterone-3-BSA Antigen FV505

**Product Name** Testosterone-3-BSA Antigen

**Catalog Number** FV505

**Tested Applications** LF

**Brand** EastCoast Bio

**Form/Appearance** Liquid; Antigen-BSA conjugate

**Storage** Short Term ( $\leq 2$  weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Buffer** 10 mM Phosphate, pH 7.0

**Preservatives** 0.05% Sodium Azide