

## Product Datasheet

### Rat TGF beta 1 ELISA Kit 670.020

<b>Product Name</b>	Rat TGF beta 1 ELISA Kit
<b>Catalog Number</b>	670.020
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	62.5 à 2000 pg/ml
<b>Sensitivity</b>	48 pg/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	TGFb1, TGF-beta1
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other rat cytokines
<b>Specificity</b>	Recognizes both natural and recombinant rat TGF-b1
<b>Sample Types</b>	Serum, Plasma, Cell culture supernatant