

---

## Product Datasheet

### Human Serum Albumin Antigen, Native Reagent Grade GF800

<b>Product Name</b>	Human Serum Albumin Antigen, Native Reagent Grade
<b>Catalog Number</b>	GF800
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	≤-20°C
<b>Note</b>	Testing: Shown to be non reactive for HBsAg, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA approved tests.
<b>Purity</b>	Reagent Grade
<b>Source</b>	Human Plasma
<b>Specificity</b>	For further testing results, please inquire
<b>Preservatives</b>	None