
Product Datasheet

Human Serum Albumin Antigen, Native Reagent Grade GF800

| | |
|------------------------|---|
| Product Name | Human Serum Albumin Antigen, Native Reagent Grade |
| Catalog Number | GF800 |
| Brand | EastCoast Bio |
| Form/Appearance | Lyophilized |
| Storage | ≤-20°C |
| Note | Testing: Shown to be non reactive for HBsAg, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA approved tests. |
| Purity | Reagent Grade |
| Source | Human Plasma |
| Specificity | For further testing results, please inquire |
| Preservatives | None |