

Product Datasheet

Human Serum Albumin Antigen, Native Reagent Grade GF800

Product Name	Human Serum Albumin Antigen, Native Reagent Grade
Catalog Number	GF800
Brand	EastCoast Bio
Form/Appearance	Lyophilized
Storage	≤-20°C
Note	Testing: Shown to be non reactive for HBsAg, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA approved tests.
Purity	Reagent Grade
Source	Human Plasma
Specificity	For further testing results, please inquire
Preservatives	None