

Product Datasheet

Benzodiazepine (Oxazepam) Antibody GN100

Product Name	Benzodiazepine (Oxazepam) Antibody
Catalog Number	GN100
Species/Host	Sheep
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, Ig Fraction
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Cross Reactivity	Benzodiazepines including Flunitrazepam, Clonazepam, 7-amino Clonazepam, Lorazepam, and Bromazepam
Specificity	Oxazepam and related compounds
Matched Pairs	GN508 (capture)
Preservatives	0.05% Sodium Azide