

## Product Datasheet

### Amphetamine Antibody GN300

**Product Name** Amphetamine Antibody

**Catalog Number** GN300

**Species/Host** Mouse

**Tested Applications** LF

**Brand** EastCoast Bio

**Form/Appearance** Purified Monoclonal Antibody

**Storage** Short Term ( $\leq 2$  weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Isotype** IgG2b, kappa

**Clonality** Monoclonal

**Buffer** 10 mM Phosphate Buffered Saline, pH 7.2

**Specificity** Amphetamine

**Preservatives** 0.05% Sodium Azide