

## Product Datasheet

### Herpes Simplex Virus 1 (HSV 1-ICP0) Antibody H1A027

<b>Product Name</b>	Herpes Simplex Virus 1 (HSV 1-ICP0) Antibody
<b>Catalog Number</b>	H1A027
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, WB, IFA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1, kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Specificity</b>	ICP0