

Product Datasheet

Herpes Simplex Virus 2 (HSV 2-ICP8) Antibody H2A024

Product Name	Herpes Simplex Virus 2 (HSV 2-ICP8) Antibody
Catalog Number	H2A024
Species/Host	Mouse
Tested Applications	ELISA, WB, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	ICP8