

Product Datasheet

Clonazepam Antibody HM003

Product Name	Clonazepam Antibody
Catalog Number	HM003
Species/Host	Mouse
Tested Applications	LF
Immunogen	Clonazepam-KLH
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Clonazepam
Matched Pairs	LA006 (capture)
Preservatives	0.1% Sodium Azide