

Product Datasheet

Fentanyl Antibody HM004

Product Name Fentanyl Antibody

Catalog Number HM004

Species/Host Mouse

Tested Applications LF

Brand EastCoast Bio

Form/Appearance Protein A Purified Monoclonal Antibody

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Clonality Monoclonal

Buffer Phosphate Buffered Saline, pH 7.4

Cross Reactivity Norfentanyl 5%, Carfentanil oxalate 0.1%

Specificity Fentanyl

Matched Pairs LA008 (capture)

Preservatives 0.05% Sodium Azide