

Product Datasheet

Ketamine Antibody HM006

Product Name Ketamine Antibody

Catalog Number HM006

Species/Host Mouse

Tested Applications LF

Immunogen Ketamine-BSA

Brand EastCoast Bio

Form/Appearance DEAE Purified Monoclonal Antibody

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG2a

Clonality Monoclonal

Purity >90%

Buffer 10 mM Potassium Phosphate, 0.85% NaCl, pH 7.2

Cross Reactivity Ketamine 100% NorKetamine: 60%

Specificity Ketamine

Matched Pairs LA010 (capture)

Preservatives 0.05% Sodium Azide