

Product Datasheet

Methamphetamine Antibody HM008

Product Name Methamphetamine Antibody

Catalog Number HM008

Species/Host Mouse

Tested Applications LF

Brand EastCoast Bio

Form/Appearance Protein A Purified Monoclonal Antibody

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG

Clonality Monoclonal

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Specificity Methamphetamine

Matched Pairs LA013 or LA014 (capture)

Preservatives 0.1% Sodium Azide