

Product Datasheet

Human IL-1 beta ELISA Kit 850.006

Product Name	Human IL-1 beta ELISA Kit
Catalog Number	850.006
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	15.6-500 pg/ml
Sensitivity	7 pg/ml
Tested Applications	ELISA
Alternative Names	IL-1b, IL-1beta
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural (pro and mature) and recombinant human IL-1b
Sample Types	Serum, Plasma, Cell culture supernatant