

Product Datasheet

Streptococcus pneumoniae Antibody HM047

Product Name	Streptococcus pneumoniae Antibody
Catalog Number	HM047
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG3
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	C-polysaccharide (phosphorylcholine epitope)
Matched Pairs	HM048 (capture)
Preservatives	0.1% Sodium Azide