

Product Datasheet

Rotavirus Antibody HM053

Product Name	Rotavirus Antibody
Catalog Number	HM053
Species/Host	Mouse
Tested Applications	LF
Immunogen	Purified human rotavirus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Rotavirus
Matched Pairs	HM053 (self-pairs) or HM052 (detection)
Preservatives	0.02% Sodium Azide