

Product Datasheet

Adenovirus Antibody HM057

Product Name	Adenovirus Antibody
Catalog Number	HM057
Species/Host	Mouse
Tested Applications	LF
Immunogen	Purified human adenovirus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	Phosphate Buffered Saline, pH 7.0
Specificity	Adenovirus
Matched Pairs	HM056 (capture)
Preservatives	0.02% Sodium Azide