

Product Datasheet

Hemoglobin (Hb) Antibody HM061

Product Name	Hemoglobin (Hb) Antibody
Catalog Number	HM061
Species/Host	Mouse
Tested Applications	LF
Immunogen	Purified Human hemoglobin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Hemoglobin
Matched Pairs	HM754 or HM753 (capture)
Preservatives	0.09% Sodium Azide