

Product Datasheet

Procalcitonin (PCT-N terminal) Antibody HM078

Product Name	Procalcitonin (PCT-N terminal) Antibody
Catalog Number	HM078
Species/Host	Mouse
Tested Applications	LF
Immunogen	Purified recombinant human Procalcitonin
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (\leq 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	PCT N-terminal epitope
Matched Pairs	HM081 or HM220 (detection)
Preservatives	0.09% Sodium Azide