

Product Datasheet

Human IL-7 ELISA Kit 850.040

Product Name	Human IL-7 ELISA Kit
Catalog Number	850.040
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	6.25-200 pg/ml
Sensitivity	3 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-7
Sample Types	Serum, Plasma, Cell culture supernatant