

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody HM084

Product Name	Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody
Catalog Number	HM084
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Native HIV antigen
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	HIV-1 p24
Matched Pairs	HM083 or JP040