

Product Datasheet

Malaria (pf-HRP-2) Antibody HM086

| | |
|----------------------------|---|
| Product Name | Malaria (pf-HRP-2) Antibody |
| Catalog Number | HM086 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Recombinant Malaria HRP-2 antigen |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | -20°C. Avoid repeated freezing and thawing. |
| Clonality | Monoclonal |
| Purity | ≥90% |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.4 |
| Specificity | Plasmodium falciparum (malaria) HRP-2 |
| Matched Pairs | HM087 (detection) |
| Preservatives | 0.1% Proclin 300 |