

## Product Datasheet

### Malaria (pf-HRP-2) Antibody HM087

<b>Product Name</b>	Malaria (pf-HRP-2) Antibody
<b>Catalog Number</b>	HM087
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, FC, Lateral Flow
<b>Immunogen</b>	Recombinant Malaria HRP-2 antigen
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>90%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Specificity</b>	Plasmodium falciparum (malaria) HRP-2
<b>Matched Pairs</b>	HM086 (capture)
<b>Preservatives</b>	0.1% Proclin 300