

Product Datasheet

Malaria (pf-pLDH) Antibody HM090

Product Name	Malaria (pf-pLDH) Antibody
Catalog Number	HM090
Species/Host	Mouse
Tested Applications	LF
Immunogen	Recombinant Malaria pLDH antigen
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4 with 0.15M NaCl
Specificity	Plasmodium falciparum pLDH
Matched Pairs	HM088 (detection)