

Product Datasheet

Human Immunoglobulin G (IgG) Antibody HM097

Product Name	Human Immunoglobulin G (IgG) Antibody
Catalog Number	HM097
Species/Host	Mouse
Tested Applications	LF
Immunogen	Human IgG
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	-20°C. Avoid repeated freezing and thawing.
Isotype	IgM
Clonality	Monoclonal
Purity	≥90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Human IgG
Preservatives	0.1% Proclin 300