

## Product Datasheet

### Human Immunoglobulin G (IgG) Antibody HM097

<b>Product Name</b>	Human Immunoglobulin G (IgG) Antibody
<b>Catalog Number</b>	HM097
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF
<b>Immunogen</b>	Human IgG
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	-20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgM
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥90%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Specificity</b>	Human IgG
<b>Preservatives</b>	0.1% Proclin 300