

## Product Datasheet

### Human Immunoglobulin M (IgM) Antibody HM098

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Human Immunoglobulin M (IgM) Antibody       |
| <b>Catalog Number</b>      | HM098                                       |
| <b>Species/Host</b>        | Mouse                                       |
| <b>Tested Applications</b> | LF, ELISA                                   |
| <b>Immunogen</b>           | Human IgM                                   |
| <b>Brand</b>               | EastCoast Bio                               |
| <b>Form/Appearance</b>     | Purified Monoclonal Antibody                |
| <b>Storage</b>             | -20°C. Avoid repeated freezing and thawing. |
| <b>Isotype</b>             | IgG1  |
| <b>Clonality</b>           | Monoclonal                                  |
| <b>Purity</b>              | >90%  |
| <b>Buffer</b>              | 10 mM Phosphate Buffered Saline, pH 7.4     |
| <b>Specificity</b>         | Human IgM                                   |
| <b>Preservatives</b>       | 0.1% Proclin 300                            |